Bedienungsanleitung
Instruction manual
Manual de instrucciones

Gay-Lussac

Pyknometer, justiert
Specific gravity bottles, calibrated
Picnómetros, ajustados

bei / at / a 20 °C ("In")

Inhalt

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1. Product description

The specific gravity bottle is a measuring device to determine the density of liquids. It is made of borosilicate glass 3.3 with a very good chemical resistance and high resistance to thermal shocks.

The specific gravity bottle acc. to Gay-Lussac consists of a pear-shaped flask and a ground-in glass stopper with capillary. The exact capacity (adjusted to contain "ln" at 20 °C) is durably engraved into the glass. The accuracy of the capacity is ± 0.001 ml.

Nominal capacities available: 1 ml, 2 ml, 5 ml, 10 ml, 25 ml, 50 ml and 100 ml

2. Measuring principle

Adjusted pycnometers according to Gay-Lussac are used to provide precisely reproducible volumes of liquids. By weighing the empty and the filled pycnometer, the density of the liquid can be calculated.

3. Handling

The pycnometer must be completely clean and dry before use. The empty mass must be determined if unknown.

Since the density is temperature-dependent, the filling and weighing must be carried out in a temperature-controlled manner, for example by temperature compensation in a water bath at the desired temperature.

The pycnometer is filled to top. When inserting the glass stopper, excess liquid penetrates through the capillary of the stopper. After filling, the pycnometer has to be dried from outside. It must be ensured that the vessel and the capillary in the glass stopper are filled completely and bubble-free. If this is the case, the filled pycnometer can be weighed.

4. Determination of density

To determine the density, the pycnometer first is weighed when empty and after that when filled and heated to 20 °C.

The density can now be calculated using the following formula:

$$\rho_{Fl} = \frac{m_2 - m_0}{m_1 - m_0} \rho_w$$

with

- m_0 being the mass of the empty pycnometer
- m_1 being the mass of the pycnometer filled with water
- m_2 being the mass of the pycnometer filled with the liquid to be examined
- ρ_w being the density of water

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